Coldex® Extra

INSTRUCTIONS FOR USE (EN)

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Coldex® Extra

Instructions for use – ENGLISH

Contents

1. INTENDED USE
2. COMPOSITION
3. INDICATIONS FOR USE
4. CONTRAINDICATIONS
5. WARNINGS
6. PRECAUTIONS
7. PRODUCT DISPOSAL
8. PRINCIPLES OF OPERATION

1. INTENDED USE

Coldex® Extra wound dressing is a moist, sterile, hydrophilic polyvinyl alcohol (PVA) foam, intended for the management of wounds via the application of Negative Pressure Wound Therapy (NPWT) with negative pressure levels between 40 mmHg and 150 mmHg. In conjunction with negative pressure wound therapy, Coldex® Extra wound dressing ensures reliable removal of wound exudate. In doing so, the applied negative pressure is evenly distributed on the wound base.

An intact sealing system guarantees a reliable airtight closure of the wound and thus an optimal use of negative pressure and exudate transport.

2. COMPOSITION

Polyvinyl alcohol

3. INDICATIONS FOR USE

Coldex® Extra can be used in skin defects, where primary closure is not possible or not desired, by reducing edema, supporting the formation of granulation tissue, promoting perfusion and by removing exudate and infectious material. The dense structure of the foam helps to restrict ingrowth of granulation tissue. Due to the higher density of Coldex® Extra a reduced level of negative pressure may be delivered through long sections of the foam.

Examples of appropriate wound types include:

- Acute
- Chronic
- Traumatic
- Sub-acute and dehisced wounds
- Ulcers such as caused by pressure (PU), diabetes (DFU) and perfusion disorders (VLU)
- Partial-thickness burns
- Flaps and grafts
- Explored undermined and tunneling wounds

Coldex® Extra should only be used by medical and nursing staff.

4. CONTRAINDICATIONS

Use of Coldex® Extra is contraindicated in presence of:

- Untreated osteomyelitis
- Exposed arteries, veins, organs or nerves

- Necrotic tissue with eschar present
- Malignancy of the wound (with exception of palliative care to enhance quality of life)
- Non-enteric and unexplored fistulas
- Exposed anastomotic sites

5. WARNINGS

1) If unexpected **blood loss** is detected, immediately discontinue therapy and take appropriate measures to control the bleeding.

Some patients (listed below) are at high risk of bleeding complications, which if uncontrolled, could be potentially fatal:

- Patients suffering from difficult hemostasis
- Patients with weakened/fragile blood vessels/organs around the wound
- Patients receiving anticoagulant therapy or platelet aggregation inhibitors
- Patients with intrinsic coagulation problems such as low platelet counts
- Patients who do not have adequate tissue coverage over vascular structures

Patients with any of the above bleeding complications must be treated and monitored in an appropriate care setting.

2) Ensure that all exposed tendons, ligaments and bones in/around the wound have adequate coverage. **Sharp edges** such as exposed bone fragments must be removed or covered as they can damage blood vessels and organs.

3) Carefully monitor the patient and dressing frequently for any signs of bleeding or signs of local or systemic infection or loss of NPWT.

4) Monitor the wound on a regular basis for signs of **infection** and make sure that all wound filler pieces are removed with each dressing change to reduce the risk of infection.

5) Monitor for **unexpected fluids** appearing in the dressing, tubing or canister such as fecal, urinary or cerebral spinal fluid; if this occurs stop therapy.

6) In general, the **weight, size and age as well as the general condition** of the patient should be taken into consideration when carrying out treatment with negative pressure wound therapy.

7) More frequent device and wound dressing monitoring should be conducted for patients who are or may be:

- Suffering from infected blood vessels
- Untreated for malnutrition
- Noncompliant or combative
- Suffering from wounds in close proximity to blood vessels or delicate fascia

8) If a patient with **spinal cord injury** experiences autonomic dysreflexia, discontinue use of NPWT.

9) Do not place NPWT in proximity of the vagus nerve to minimize risk of **bradycardia**.

10) Coldex® Extra must not be packed too tightly or forced into any wound area.

11) Avoid use of Coldex® Extra **circumferentially** except in cases of edema or heavily exuding extremities, where this technique may be necessary to maintain a seal. Consider using multiple drapes to minimize risk of decreased distal pulses and discontinue therapy if changes in circulation are detected.

12) Do not place foam into blind or unexplored tunnels. If a tunnel of known depth is present, cut a wedge shaped piece of foam, which is longer than the tunnel to ensure direct contact with the foam in the primary wound cavity.

13) If **defibrillation** or **MRI scan** is required, the foam dressing can be left in place if the location of the dressing does not interfere with the treatment. Other devices for the therapy need to be disconnected. When **hyperbaric oxygen units** are used, the dressing must be replaced with a dressing that is compatible with the hyperbaric oxygen environment. Never stop therapy for more than 2 hours.

6. PRECAUTIONS

- 1) Do not use Coldex® Extra if the package is breached or damaged.
- 2) Coldex® Extra is supplied sterile; do not re-sterilize.
- 3) Do not apply alcoholic solutions on Coldex® Extra as this compromises the integrity, composition and functionality of the foam.
- 4) Coldex® Extra is for single use only; reuse of the device can lead to cross contamination.
- 5) Coldex® Extra is supplied moist and should be used as soon as possible after opening.
- 6) Coldex® Extra does not contain antibiotics. It has no bacteriostatic effect on pre-existing infections, nor does it prevent the development of new infections.
- 7) Do not cut Coldex® Extra directly over the wound cavity to avoid foam fragments from falling into the wound. Rub edges of the foam away from the open wound.
- 8) Always count the total number of dressings used to pack a wound. Document the number in the patient's records to ensure all the dressings are accounted for upon removal.
- 9) Do not use Coldex® Extra after the expiry date.
- 10) Keep Coldex® Extra away from sunlight and heat; store Coldex® Extra between 5°C and 35°C. Keep Coldex® Extra dry during storage.

7. PRODUCT DISPOSAL

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of Coldex® Extra should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

8. PRINCIPLES OF OPERATION


Use clean or aseptic techniques for application according to your institutional protocol. Thorough wound assessment and cleansing should occur with each dressing change. Assess the patient's wound characteristics including size and shape, wound exudate levels, any undermining/tunneling, the presence of any veins, arteries, organs, tendons, ligaments or nerves and osteomyelitis.

Asses, cleanse and debride the wound

Assess the wound dimensions and characteristics. Debride any necrotic or eschar tissue. Cleanse the wound bed according to your institutional protocol. Cleansing of the wound must be performed before applying the dressing for the first time and after every dressing change.

Dress the wound

Ensure that all exposed bone, tendons and ligaments in/around the wound are completely covered and protected by an intermediate, non-stick wound cushion (limiting membrane). Cut Coldex® Extra with sterile instruments to fit the size and shape of the wound. Do not cut Coldex® Extra over the wound, as fragments may fall into the wound bed. Always rub Coldex® Extra away from the wound bed to remove loose edges.

Gently place Coldex® Extra into the wound. Make sure Coldex® Extra does not cover the skin around the wound or growing epithelium at the wound edge.


Count number of pieces used in the wound. If multiple pieces are needed, count and record how many pieces are used to ensure all pieces are removed with a dressing change to minimize the risk of retention and possible infection.

Avoid over packing or forcing the dressing into the wound.

If using in an explored tunneling wound, ensure Coldex® Extra is longer than the tunnel and wedge shaped to assist easy removal.

Seal the wound

An intact sealing system guarantees a reliable airtight closure of the wound and thus an optimal use of negative pressure and exudate transport.

Always read the user manual(s) of the relevant device(s) when setting up a NPWT system to ensure a safe course of therapy.

Dressing changes

Each time the dressing is changed the wound must be thoroughly cleansed in accordance with the physician's instructions and applicable hygiene standards. If the patient complains of pain while the dressing is being changed, premedication or the use of a non-adhering wound contact layer (e.g. silicone layer) or the administration of a local anesthetic should be considered.

Coldex® Extra should be changed every 48 to 72 hours after the initial start of therapy. If no leakage is present and the patient is in a comfortable state, dressing changes should be performed at least 3 times per week. Frequent inspection is recommended and more frequent dressing changes may be required.

If there is evidence of an infection of the wound, Coldex® Extra must be changed regularly, at least once every 24 hours until the signs of the infection have decreased significantly or disappeared altogether.

If dressing adheres to the wound, apply saline or sterile water to the wound dressing. Once the dressing is fully hydrated the foam can be removed.

Apply and remove dressing carefully from sensitive or fragile skin to avoid skin stripping, especially after frequent dressing changes. Use of a skin sealant may assist with protection of periwound skin.

Count number of pieces taken out of the wound. Check patient's record to see how many pieces were used in order to ensure all pieces are removed to minimize the risk of retention and possible infection.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, who will take care for reporting the applicable Competent Authority.

 Mondomed NV
Middenweg 20
BE-3930 Hamont-Achel

 1639

 +32 11 440 990
Email: info@mondomed.be
<http://www.mondomed.be>

	European conformity		Do not use if package is damaged
	Medical device		Distributor
	Manufacturer		Consult instructions for use
	Date of manufacture		Caution
	Use-by date		Keep away from sunlight
	Batch code		Keep dry
	Catalogue number		Temperature limit
	Sterilized using irradiation		Green dot
	Single sterile barrier system		Dispose of packaging in an environmentally responsible way
	Unique device identification		Composite of Paper and plastic/aluminum
	Do not re-use		Non-corrugated fiberboard (paperboard)
	Do not re-sterilize		Corrugated fiberboard (cardboard)

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