Fractomed[®] Lite

instructions for use – ENGLISH

Contents

1. INTENDED USE
2. COMPOSITION
3. INDICATIONS FOR USE
4. CONTRAINDICATIONS
5. PRECAUTIONS
6. PRODUCT DISPOSAL
7. PRINCIPLES OF OPERATION

1. INTENDED USE

Fractomed[®] Lite is a rapid response ready-to-use splinting material for fast and accurate first aid and emergency interventions and as immobilization after postoperative treatment. Its microporosity ensures that sufficient perspiration is possible and prevents decubitus, swelling, skin inflammation, etc. Fractomed[®] Lite is radiolucent which means no need to remove the device to make an x-ray. Fractomed[®] Lite is available in different sizes and forms and in left and right anatomical forms.

2. COMPOSITION

Polyvinyl alcohol (PVA) foam

3. INDICATIONS FOR USE

Non exhaustive list

- Fractures
- Dislocations
- Sprains
- Muscle tendon contractures
- Postoperative treatment
- First aid

Fractomed[®] Lite should only be used by medical professionals.

4. CONTRAINDICATIONS

Fractomed[®] may not be used directly on open wounds.

5. PRECAUTIONS

- 1) Do not apply alcoholic solutions on the device as this compromises the integrity, composition

and functionality of the foam.

- 2) Due to the properties of the material and the properties of the bandage, Fractomed[®] Lite can easily be adjusted if it is too tight or too loose.
- 3) Do not use the device after the expiry date.
- 4) Keep the device away from sunlight and heat. Keep Fractomed[®] Lite dry during storage.

6. PRODUCT DISPOSAL

In order to minimize the risk of potential infection hazards, or environmental pollution disposal of Fractomed[®] Lite should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

7. PRINCIPLES OF OPERATION

Apply a (thin) stockinet on the limb to be treated. If necessary additional padding material can be used. Take the splint out of the package and position it on the limb. Check that the dimensions and shapes are correct and where necessary make adjustments. The material can be cut with a plaster scissors. Fix the splint in its position with the semi-elastic bandage.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, who will take care for reporting the applicable Competent Authority.

 Mondomed NV
Middenweg 20
BE-3930 Hamont-Achel
Phone: +32 11 440 990
Email: info@mondomed.be

<http://www.mondomed.be>

	European conformity
	Medical device
	Manufacturer
	Use-by date
	Batch code
	Catalogue number
	Unique Device Identifier
	Single patient - multiple use
	Keep away from sunlight
	Keep dry
	Dispose of packaging in an environmentally responsible way
	Green dot

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