## **RetractorCover**

INSTRUCTIONS FOR USE (EN)

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## RetractorCover

Instructions for use – ENGLISH

### Contents

1. INTENDED USE
2. COMPOSITION
3. INDICATIONS FOR USE
4. CONTRAINDICATIONS
5. WARNINGS
6. PRECAUTIONS
7. PRODUCT DISPOSAL
8. PRINCIPLES OF OPERATION

#### 1. INTENDED USE

If during surgery, point and blade of the retractor are not protected damage of tissues and organs will be possible. RetractorCover is an appropriate solution to this kind of inconvenience as it is designed to offer protection to all tissues and organs. In addition, reflection of light is reduced and the stability of the retractor will be enhanced because it does not slip away after insertion.

RetractorCover has a porous structure with good absorbent properties.

#### 2. COMPOSITION

Polyvinyl alcohol (PVA) foam  
X-ray detectable barium thread

#### 3. INDICATIONS FOR USE

Surgical procedures where a retractor is used.

RetractorCover should only be used by medical and nursing staff.

#### 4. CONTRAINDICATIONS

No known contraindications.

#### 5. WARNINGS

- 1) Before use, the RetractorCover must be moistened by immersing it in sterile saline for approximately 20 seconds.
- 2) RetractorCover must be kept moist during entire duration of use. Prevent RetractorCover from drying out during the procedure.
- 3) RetractorCover must be removed after surgery. RetractorCover is equipped with an X-ray detectable tread and can be visualized with X-ray if necessary.

#### 6. PRECAUTIONS

- 1) Do not use the RetractorCovers if package is breached or damaged.
- 2) RetractorCovers are supplied sterile; do not re-sterilize.
- 3) Do not apply alcoholic solutions on the RetractorCover as this can compromise the integrity, composition and functionality of the foam.
- 4) RetractorCover is for single use only. Reuse of the device can lead to several hazards:

- Cross contamination
- Infection risks
- Device failure caused by initial use resulting in patient injury

5) Do not use RetractorCovers after the expiry date.

6) Keep RetractorCover away from sunlight and heat. Keep dry during storage and store RetractorCovers between 5°C and 35°C.

## 7. PRODUCT DISPOSAL

In order to minimize the risk of potential infection hazards, or environmental pollution, RetractorCover should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

## 8. PRINCIPLES OF OPERATION

Use clean or aseptic techniques for application, according to your institutional protocol.

Before use, the RetractorCover must be moistened by immersing it in sterile saline for approximately 20 seconds.

Once the RetractorCover is saturated, squeeze out the excess saline. The RetractorCover is now ready for use and can be slipped over the retractor.

The RetractorCover must kept moist while in use and need to be removed at end of procedure.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, who will take care for reporting the applicable Competent Authority.

 Mondomed NV  
Middenweg 20  
BE-3930 Hamont-Achel

 1639

 +32 11 440 990  
Email: [info@mondomed.be](mailto:info@mondomed.be)  
<http://www.mondomed.be>

	European conformity		Do not resterilize
	Medical device		Do not use if package is damaged
	Manufacturer		Do not re-use
	Date of manufacture		Consult instructions for use
	Use-by date		Keep away from sunlight
	Batch code		Keep dry
	Catalogue number		Temperature limit

Sterilized using irradiation

Green dot

Double sterile barrier system

Dispose of packaging in an environmentally responsible way

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