Mondomed Retractorpad®

INSTRUCTIONS FOR USE (EN)

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BRUGSANVISNING (DA)

BRUKSANVISNING (NO)

BRUKSANVISNING (SV)

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Instructions for use – ENGLISH

Contents

1. INTENDED USE
2. COMPOSITION
3. INDICATIONS FOR USE
4. CONTRAINDICATIONS
5. WARNINGS
6. PRECAUTIONS
7. PRODUCT DISPOSAL
8. PRINCIPLES OF OPERATION

1. INTENDED USE

Retractorpad is intended to be used in abdominal and pelvic surgery to keep organs and tissues, in particular the intestines out of the field of surgery.

Retractorpad is slightly moist, flexible and has a soft texture preventing damage or contact trauma to the organs and tissues, in particular the intestines and peritoneum. As the Retractorpad is moist, dehydration of tissues is avoided. Because of its spongy, porous structure blood and other body fluids will be absorbed.

Retractorpad has its support on the paracolic gutters thanks to its U-shaped structure, while the arche of the device prevents the underlying blood vessels from being compressed.

2. COMPOSITION

Polyvinyl alcohol (PVA) foam

X-ray detectable barium thread

3. INDICATIONS FOR USE

Retractorpad is used in abdominal and pelvic surgical procedures.

Retractorpad should only be used by medical and nursing staff.

4. CONTRAINDICATIONS

No known contraindications.

5. WARNINGS

1) Retractorpad must be kept moist during entire duration of use.

2) Make sure the large blood vessels next to the vertebral column are not compressed.

3) After surgery, Retractorpad should be removed. In case Retractorpad can't be found it can be detected by the X-ray visible thread.

6. PRECAUTIONS

1) Do not use Retractorpad if package is breached or damaged.

2) Retractorpad is supplied sterile; do not re-sterilize.

3) Do not apply alcoholic solutions on the pads as this compromises the integrity, composition and functionality of the foam.

4) Mondomed Retractorpad® is for single use only; reuse of the device can lead to several hazards:

- Cross contamination

- Infection risks
- Deterioration and alteration of the characteristics of the device due to the first use with consequent risk of injury to the patient

5) Do not use the Retractorpad after the expiry date.

6) Keep the Retractorpad away from sunlight and heat. Keep the pads dry during storage and store them at temperatures between 5°C and 35°C.

7. PRODUCT DISPOSAL

In order to minimize the risk of potential infection hazards, or environmental pollution Retractorpad should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

8. PRINCIPLES OF OPERATION

Introduce the double-folded Retractorpad in the surgical area while using one hand to hold the Retractorpad and the other hand to hold the intestines out of the field of surgery.

If necessary, Retractorpad can be cut to optimal size. Do not cut Retractorpad over the surgical area, as fragments may fall into the wound bed.

The double-folded Retractorpad is then unfolded in order to install the notch over the vertebral column, while the two paw-shaped configurations of the Retractorpad can be placed on the left and the right side of the vertebral column.

In order to verify whether there is any compression of the large blood vessels next to the vertebral column, there must be ample space to place one finger between the U-shaped notch of the lower edge of the Retractorpad and the vertebral column.

It is sometimes necessary to use a speculum in order to keep the Retractorpad in its place. The Retractorpad can be used together with any type of abdominal spreader.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, who will take care for reporting the applicable Competent Authority.

 Mondomed NV
Middenweg 20
BE-3930 Hamont-Achel

 1639

 +32 11 440 990
Email: info@mondomed.be
<http://www.mondomed.be>

	European conformity		Do not resterilize
	Medical device		Do not use if package is damaged
	Manufacturer		Consult instructions for use
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Temperature limit

Catalogue number

Do not re-use

Sterilized using irradiation

Green dot

Double sterile barrier system

Dispose of packaging in an environmentally responsible way

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